

Product Datasheet

ST2 (IL1RL1) antibody (orb1327999)

Description	Mouse monoclonal antibody to ST2 (IL1RL1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full-length human recombinant protein of human ST2(NP_057316) produced in HEK293E
Target	IL1RL1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	63.36 kDa
Uniprot ID	Q01638
NCBI	NM_016232
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.