

Product Datasheet

Serotonin N acetyltransferase (AANAT) antibody (orb1328015)

Description	Mouse monoclonal antibody to Serotonin N acetyltransferase (AANAT)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-72 of human AANAT (NP_001079) produced in E.coli.
Target	AANAT
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	23.2 kDa
Uniprot ID	Q16613
NCBI	NM_001088
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

