


Product Datasheet

SHARP2 (BHLHE40) antibody (orb1328038)

Description	Mouse monoclonal antibody to SHARP2 (BHLHE40)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-313 of human BHLHE40 (NP_003661) produced in E.coli.
Target	BHLHE40
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	45.3 kDa
Uniprot ID	O14503
NCBI	NM_003670
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...