

Product Datasheet

NeuN (RBFOX3) antibody (orb1328159)

Description	Mouse monoclonal antibody to NeuN (RBFOX3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-97 of human RBFOX3(NP_001076044) produced in E.coli.
Target	RBFOX3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	33.7 kDa
Uniprot ID	A6NFN3
NCBI	NM_001082575
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.

