

Product Datasheet

Ceramide synthase 2 (CERS2) antibody (orb1328182)

Description Mouse monoclonal antibody to Ceramide synthase 2 (CERS2)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Human recombinant protein fragment corresponding to amino acids 62-136 of human CERS2(NP_859530) produced in E.coli.

Target CERS2

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2a

Clonality Monoclonal


MW 44.7 kDa

Uniprot ID [Q96G23](#)

NCBI [NM_181746](#)

Dilution Range WB 1:500~2000, IHC 1:1000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...