


Product Datasheet

MAG1 (GPAT3) antibody (orb1328266)

Description	Mouse monoclonal antibody to MAG1 (GPAT3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 41-131 of human AGPAT9 (NP_116106) produced in E.coli.
Target	GPAT3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	48.5 kDa
Uniprot ID	Q53EU6
NCBI	NM_032717
Dilution Range	WB 1:500, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...