


Product Datasheet

GNMB antibody (orb1328345)

Description	Mouse monoclonal antibody to GNMB
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 508-572 of human GNMB (NP_001005340) produced in E.coli.
Target	GNMB
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	61.5 kDa
Uniprot ID	Q14956
NCBI	NM_001005340
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...