


Product Datasheet

PASD5 (NPAS1) antibody (orb1328362)

Description	Mouse monoclonal antibody to PASD5 (NPAS1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 319-522 of human NPAS1 (NP_002508) produced in E.coli.
Target	NPAS1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	62.5 kDa
Uniprot ID	Q99742
NCBI	NM_002517
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...