


Product Datasheet

COL5A2 antibody (orb1328365)

Description	Mouse monoclonal antibody to COL5A2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1200-1499 of human COL5A2(NP_000384) produced in E.coli.
Target	COL5A2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	144.7 kDa
Uniprot ID	P05997
NCBI	NM_000393
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...