



Product Datasheet

MCP5 antibody (orb13284)

Description	MCP5 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from mouse CCL12 (23-98/104 aa)
Target	CCL12
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	9 kDa
Uniprot ID	Q62401
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Flow cytometric analysis of Mouse spleen cells.