

Product Datasheet

TRMT12 antibody (orb1328400)

Description	Mouse monoclonal antibody to TRMT12
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TRMT12 (NP_060426) produced in E.coli.
Target	TRMT12
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	50.1 kDa
Uniprot ID	Q53H54
NCBI	NM_017956
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

