


Product Datasheet

BCL7B antibody (orb1328449)

Description	Mouse monoclonal antibody to BCL7B
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human BCL7B (NP_001698) produced in E.coli.
Target	BCL7B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	22 kDa
Uniprot ID	Q9BQE9
NCBI	NM_001707
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...