

Product Datasheet

C14orf50 (PPP1R36) antibody (orb1328465)

Description Mouse monoclonal antibody to C14orf50 (PPP1R36)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human PPP1R36(NP_758953) produced in E.coli.

Target PPP1R36

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal

MW 49.3 kDa

Uniprot ID [Q96LQ0](#)

NCBI [NM_172365](#)

Dilution Range WB 1:500~2000

Expiration Date 12 months from date of receipt.

