

Product Datasheet

n-Myc (MYCN) antibody (orb1328548)

Description	Mouse monoclonal antibody to n-Myc (MYCN)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 229-464 of human MYCN(NP_005369) produced in E.coli.
Target	MYCN
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	49.4 kDa
Uniprot ID	P04198
NCBI	NM_005378
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

