

Product Datasheet

C1orf83 (TCEANC2) antibody (orb1328557)

Description Mouse monoclonal antibody to C1orf83 (TCEANC2)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human TCEANC2(NP_694580) produced in E.coli.

Target TCEANC2

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 24 kDa

Uniprot ID [Q96MN5](#)

NCBI [NM_153035](#)

Dilution Range WB 1:500~2000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...