


Product Datasheet

CEP63 antibody (orb1328592)

Description	Mouse monoclonal antibody to CEP63
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human CEP63 (NP_001035859) produced in E.coli.
Target	CEP63
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	62.8 kDa
Uniprot ID	Q96MT8
NCBI	NM_001042400
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...