

Product Datasheet

PATZ1 antibody (orb1328602)

Description	Mouse monoclonal antibody to PATZ1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PATZ1(NP_114440) produced in E.coli.
Target	PATZ1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	57.4 kDa
Uniprot ID	Q9HBE1
NCBI	NM_032051
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.

