


Product Datasheet

UEVLD antibody (orb1328635)

Description	Mouse monoclonal antibody to UEVLD
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 131-379 of human UEVLD(NP_060784) produced in E.coli.
Target	UEVLD
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	42.1 kDa
Uniprot ID	Q8IX04
NCBI	NM_018314
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...