


Product Datasheet

ASS1 antibody (orb1328696)

Description	Mouse monoclonal antibody to ASS1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to residues near C-terminus of human ASS1 (NP_000041).
Target	ASS1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	46.3 kDa
Uniprot ID	P00966
NCBI	NM_000050
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...