

Product Datasheet

C7orf16 (PPP1R17) antibody (orb1328786)

Description Mouse monoclonal antibody to C7orf16 (PPP1R17)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human C7ORF16 (NP_006649) produced in E.coli.

Target PPP1R17

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal

Uniprot ID [O96001](#)

NCBI [NM_006658](#)

Dilution Range WB 1:2000

Expiration Date 12 months from date of receipt.

