


Product Datasheet

PDRG1 antibody (orb1328837)

Description	Mouse monoclonal antibody to PDRG1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 33-117 of human PDRG1(NP_110442) produced in E.coli.
Target	PDRG1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	15.3 kDa
Uniprot ID	Q9NUG6
NCBI	NM_030815
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...