

Product Datasheet

Homeo box C10 (HOXC10) antibody (orb1328846)

Description	Mouse monoclonal antibody to Homeo box C10 (HOXC10)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-262 of human HOXC10(NP_059105) produced in E.coli.
Target	HOXC10
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	37.9 kDa
Uniprot ID	Q9NYD6
NCBI	NM_017409
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

