


Product Datasheet

MAD4 (MXD4) antibody (orb1328855)

Description	Mouse monoclonal antibody to MAD4 (MXD4)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human MXD4(NP_006445) produced in E.coli.
Target	MXD4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	23.3 kDa
Uniprot ID	Q14582
NCBI	NM_006454
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...