


Product Datasheet

PSF3 (GINS3) antibody (orb1328899)

Description	Mouse monoclonal antibody to PSF3 (GINS3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human GINS3(NP_073607) produced in E.coli.
Target	GINS3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	24.4 kDa
Uniprot ID	Q9BRX5
NCBI	NM_022770
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...