


Product Datasheet

SDCCAG10 (CWC27) antibody (orb1328972)

Description	Mouse monoclonal antibody to SDCCAG10 (CWC27)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 145-472 of human CWC27(NP_005860) produced in E.coli.
Target	CWC27
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	53.7 kDa
Uniprot ID	Q6UX04
NCBI	NM_005869
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...