


Product Datasheet

CCL22 antibody (orb13290)

Description	CCL22 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr
Immunogen	KLH conjugated synthetic peptide derived from mouse C-C motif chemokine 22 (51-93/93 aa)
Target	CCL22
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	10 kDa
Uniprot ID	O88430
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat lung...