

Product Datasheet

TFIIB (GTF2B) antibody (orb1329035)

Description	Mouse monoclonal antibody to TFIIB (GTF2B)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 43-316 of human GTF2B(NP_001505) produced in E.coli.
Target	GTF2B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	34.7 kDa
Uniprot ID	Q00403
NCBI	NM_001514
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

