


Product Datasheet

Amelotin (AMTN) antibody (orb1329079)

Description	Mouse monoclonal antibody to Amelotin (AMTN)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human AMTN (NP_997722) produced in E.coli.
Target	AMTN
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	21.4 kDa
Uniprot ID	Q6UX39
NCBI	NM_212557
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...