


Product Datasheet

EXOSC1 antibody (orb1329082)

Description	Mouse monoclonal antibody to EXOSC1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human EXOSC1 (NP_057130) produced in E.coli.
Target	EXOSC1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	21.3 kDa
Uniprot ID	Q9Y3B2
NCBI	NM_016046
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...