

Product Datasheet

C11orf17 (AKIP1) antibody (orb1329098)

Description Mouse monoclonal antibody to C11orf17 (AKIP1)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC

Immunogen Full length human recombinant protein of human AKIP1 (NP_878904) produced in E.coli.

Target AKIP1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 22.9 kDa

Uniprot ID [Q9NQ31](#)

NCBI [NM_182901](#)

Dilution Range IHC 1:150

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...