


Product Datasheet

SFRS9 (SRSF9) antibody (orb1329138)

Description	Mouse monoclonal antibody to SFRS9 (SRSF9)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SRSF9 (NP_003760) produced in E.coli.
Target	SRSF9
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	25.4 kDa
Uniprot ID	Q13242
NCBI	NM_003769
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...