


Product Datasheet

Neuropeptide Y (NPY) antibody (orb1329232)

Description	Mouse monoclonal antibody to Neuropeptide Y (NPY)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NPY (NP_000896) produced in E.coli.
Target	NPY
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	8 kDa
Uniprot ID	P01303
NCBI	NM_000905
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...