


Product Datasheet

Cyclin Y (CCNY) antibody (orb1329312)

Description	Mouse monoclonal antibody to Cyclin Y (CCNY)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CCNY (NP_859049) produced in HEK293T cell.
Target	CCNY
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q8ND76
NCBI	NM_181698
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...