

Product Datasheet

GM CSF (CSF2) antibody (orb1329364)

Description	Mouse monoclonal antibody to GM CSF (CSF2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 18-144 of human CSF2(NP_000749) produced in E.coli.
Target	CSF2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	14.4 kDa
Uniprot ID	P04141
NCBI	NM_000758
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

