

Product Datasheet

Cardiac Troponin T (TNNT2) antibody (orb1329400)

Description Mouse monoclonal antibody to Cardiac Troponin T (TNNT2)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human TNNT2 (NP_000355) produced in E.coli.

Target TNNT2

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 35.4 kDa

Uniprot ID [P45379](#)

NCBI [NM_000364](#)

Dilution Range WB 1:2000, IHC 1:150

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...