

Product Datasheet

IL3RB (CSF2RB) antibody (orb1329441)

Description Mouse monoclonal antibody to IL3RB (CSF2RB)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC

Immunogen Full length human recombinant protein of human CSF2RB (NP_000386) produced in HEK293T cell.

Target CSF2RB

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 95.7 kDa

Uniprot ID [P32927](#)

NCBI [NM_000395](#)

Dilution Range IHC 1:150

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...