


Product Datasheet

TNF alpha (TNF) antibody (orb1329471)

Description	Mouse monoclonal antibody to TNF alpha (TNF)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human TNF (NP_000585) produced in HEK293T cell.
Target	TNF
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P01375
NCBI	NM_000594
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...