


Product Datasheet

TCEAL1 antibody (orb1329563)

Description	Mouse monoclonal antibody to TCEAL1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TCEAL1 (NP_004771) produced in E.coli.
Target	TCEAL1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	18.5 kDa
Uniprot ID	Q15170
NCBI	NM_004780
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...