

Product Datasheet

BMF antibody (orb1329570)

Description	Mouse monoclonal antibody to BMF
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human BMF (NP_001003940) produced in E.coli.
Target	BMF
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	20.3 kDa
Uniprot ID	Q96LC9
NCBI	NM_001003940
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

