


Product Datasheet

TRIB3 antibody (orb1329615)

Description	Mouse monoclonal antibody to TRIB3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-200 of human TRIB3(NP_066981) produced in E.coli.
Target	TRIB3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	39.4 kDa
Uniprot ID	Q96RU7
NCBI	NM_021158
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...