

Product Datasheet

FBXO8 antibody (orb1329651)

Description	Mouse monoclonal antibody to FBXO8
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-231 of human FBXO8(NP_036312) produced in E.coli.
Target	FBXO8
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	36.9 kDa
Uniprot ID	Q9NRD0
NCBI	NM_012180
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

