


Product Datasheet

KCHIP2 (KCNIP2) antibody (orb1329699)

Description	Mouse monoclonal antibody to KCHIP2 (KCNIP2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human KCNIP2 (NP_775284) produced in E.coli.
Target	KCNIP2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	28.8 kDa
Uniprot ID	Q9NS61
NCBI	NM_173192
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...