

Product Datasheet

SFRP2 antibody (orb1329736)

Description	Mouse monoclonal antibody to SFRP2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 25-295 of human SFRP2(NP_003004) produced in E.coli.
Target	SFRP2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	31.4 kDa
Uniprot ID	Q96HF1
NCBI	NM_003013
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.

