


Product Datasheet

KLF2 antibody (orb1329924)

Description	Mouse monoclonal antibody to KLF2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human KLF2 (NP_057354) produced in E.coli.
Target	KLF2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	37.2 kDa
Uniprot ID	Q9Y5W3
NCBI	NM_016270
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...