


Product Datasheet

PD1 (PDCD1) antibody (orb1330005)

Description	Mouse monoclonal antibody to PD1 (PDCD1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, NeA, WB
Immunogen	Full length human recombinant protein of human PDCD1 (NP_005009) produced in HEK293T cell.
Target	PDCD1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	29.2 kDa
Uniprot ID	Q15116
NCBI	NM_005018
Dilution Range	WB 1:2000, IHC 1:150, IF 1:900, FLOW 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...