


Product Datasheet

CCN2 antibody (orb1330011)

Description	Mouse monoclonal antibody to CCN2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 27-349 of human CTGF(NP_001892) produced in E.coli.
Target	CCN2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	37.9 kDa
Uniprot ID	P29279
NCBI	NM_001901
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...