

Product Datasheet

Desmuslin (SYNM) antibody (orb1330062)

Description	Mouse monoclonal antibody to Desmuslin (SYNM)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 329-580 of human SYNM(NP_663780) produced in E.coli.
Target	SYNM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	172.7 kDa
Uniprot ID	O15061
NCBI	NM_145728
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...