


Product Datasheet

UBE2D2 antibody (orb1330083)

Description	Mouse monoclonal antibody to UBE2D2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human UBE2D2 (NP_003330) produced in E.coli.
Target	UBE2D2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	16.6 kDa
Uniprot ID	P62837
NCBI	NM_003339
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...