

Product Datasheet

BCL2 antibody (orb1330087)

Description	Mouse monoclonal antibody to BCL2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human BCL2 (NP_000624) produced in HEK293T cell.
Target	BCL2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P10415
NCBI	NM_000633
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

