


Product Datasheet

ARID1A antibody (orb1330164)

Description	Mouse monoclonal antibody to ARID1A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1986-2285 of human ARID1A (NP_006006) produced in E.coli.
Target	ARID1A
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	241.9 kDa
Uniprot ID	O14497
NCBI	NM_006015
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...